



Instrument reprocessing in Miele dental washer-disinfectors




After treatment


- Remove cement residue directly after treatment
- Store instruments in their dry state prior to cleaning/disinfection
- Process instruments within 6 hours after use (observe manufacturer's recommendation)



- Load washer-disinfector
- Place lumened instruments on injector nozzles or adapters for hand-pieces and contra-angles
- Place hand instruments in mesh trays
- Place trays, beakers in inserts



- Make sure not to block spray arms
- Avoid positioning items in the shadow of others
- Only reprocess corrosion-proof and thermally stable instruments in a washer-disinfector
- Mouth mirrors should not come into contact with each other



Dispensing of process chemicals

- Note when replacing canisters: Never mix up siphons for detergent and neutralising agent!
- Dispensing of detergent, neutralising agent and, where applicable, surfactant is automatic
- Note canister fill level and replace canisters when prompted



Program selection

- Selection via program selectors 1-3 or menu for further programs
- Vario TD Dental = Standard program

End of program (depending on model version)

- PG 8581 (EcoDry):** Machine door opens automatically at end of program, allowing steam to escape
- PG 8591 (DryPlus):** Drying process is automatically activated and load is dried; optical signal on completion of drying

Checks

- Empty and clean filter in cabinet
- If fitted: Clean A 800 tubular filter in upper basket with A 804 brush at least once per week
- Direct coupling of upper basket to water supply
- Check cleaning performance
- Disinfection parameters
- Check non-return valve for debris according to operating instructions

Notes:

- Note detailed instructions in Operating Instructions and special notes on the reprocessing of transmission instruments